


K9 Mask® Size Chart



How to Choose the Right Size K9 Mask® for your Dog

Our K9 Mask® muzzle sizing chart will help you find the right size for your pet.

- If you do not have a flexible measuring tape you can use a shoe lace or phone charging cable to make your measurements. Then use a ruler or metal measuring tape to get the measurement.
- NOTE: If your dog is in between circumference sizes we recommend getting the larger size mask for better panting space in the mask.



K9 MASK® SIZE	DOG BREED EXAMPLES (Measure for Accuracy)	MUZZLE CIRCUMFERENCE RANGE (Most Important)	MUZZLE LENGTH RANGE
Small	Examples: Yorkshire Terrier, Shitzu, Miniature Pinscher, Ratter, Miniature Schnauzer, etc.	4.5 - 6.5 inches (11.5 - 16.5 cm)	1.6 - 2.8 inches (4.0 - 7.0 cm)
Medium	Examples: Spaniel, Schnauzer, Fox terrier, Pinscher, Corgi, etc.	7.0 - 9.0 inches (17.5 - 23.0 cm)	2.8 - 4.3 inches (7.0 - 11.0 cm)
Large	Examples: Pointer, Dalmatian, Setter, Labrador, Retriever, Siberian Husky, etc.	9.0 - 13.0 inches (23.0 - 33.0 cm)	3.7 - 4.7 inches (9.5 - 12.0 cm)
Extra Large	Examples: German Shepherd (big), Newfoundland, Mastiff, Great Dane, etc.	13.0 - 17.5 inches (33.0 - 45.0 cm)	4.5 - 5.9 inches (11.5-15.0 cm)

IMPORTANT NOTE: If your dog is between 'Muzzle Circumference' sizes choose the larger size for better panting space in the mask.





SMALL K9 MASK®

- Breed Examples: Yorkshire Terrier, Shih Tzu, Miniature Pinscher, Miniature Schnauzer. (Breeds which tend to be too small to wear this size mask: Toy Chihuahua, Brussels Griffon, Papillon.)
- Muzzle Length: 1.6-2.8 inches (4.0-7.0 cm)
- Muzzle Circumference: 4.5-6.5 inches (11.5-16.5 cm)
- ***Note:** If your dog's muzzle circumference is in between sizes we recommend you get the larger mask size for more panting space.*
- View [Small K9 Mask®](#).



MEDIUM K9 MASK®

- Breed Examples: Spaniel, Schnauzer, Fox Terrier, Pinscher.
- Muzzle Length: 2.8-4.3 inches (7-11 cm)
- Muzzle Circumference: 7.0-9.0 inches (17.5-23.0 cm)
- **Note:** *If your dog's muzzle circumference is in between sizes we recommend you get the larger mask size for more panting space.*
- View [Medium K9 Mask®](#).



LARGE K9 MASK®

- Breed Examples: German Shepherd (small), Pointer, Dalmatian, Setter, Labrador, Retriever, Husky.
- Muzzle Length: 3.7-4.7 inches (9.5-12 cm)
- Muzzle Circumference: 9.0-13.0 inches (23.0-33.0 cm)
- **Note:** *If your dog's muzzle circumference is in between sizes we recommend you get the larger mask size for more panting space.*
- View [Large K9 Mask®](#).



EXTRA LARGE K9 MASK®

- Breed Examples: German Shepherd (big), Newfoundland, Mastiff.
- Muzzle Length: 4.5-5.9 inches (11.5-15.0 cm)
- Muzzle Circumference: 13.0-17.5 inches (33.0-45.0 cm)

- ***Note:** If your dog's muzzle circumference is in between sizes we recommend you get the larger mask size for more panting space.*
- View [Extra Large K9 Mask®](#).



FLAT FACED (*Brachycephalic*) DOGS

- We **do not currently** have K9 Mask® sizes for flat faced (*Brachycephalic*) dogs like (Pug, Pekingese, Bulldog, etc.). However, we continue our research and development to find solutions for these dog breeds.